# 2x LAMP Master Mix

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Product Code: NE6041 / NE6043

## I Product Description

2x LAMP Master Mix is a reagent for Loop-mediated Isothermal Amplification (LAMP) method. This product contains isothermally amplifiable strand displacement DNA polymerase, inorganic pyrophosphatase, Mg<sup>2+</sup>, dNTPs, optimum buffer, and fluorescent dye that binds to the double-stranded DNA and allows detection of amplified DNA by fluorescence device.

### Precaution

- This product is for research use only and not for medical or other purposes.
- This product should be handled in accordance with the descriptions in the manual. We are not responsible for problems caused if this product is not handled in accordance with the manual.
- Expire date is printed on the label.
- Follow the instructions from your region to discard this product.
- This product is for fluorescence device and not for Turbidity measuring Real-time device for LAMP such as Loopamp EXIA<sup>®</sup>, LA-320C, and RT-160C (Eiken Chemical Company).
- LAMP is patented by Eiken Chemical Company. We are permitted to manufacture and sell fluorescent detecting LAMP reagents for research use by Eiken Chemical Company.

## II Contents

	100 reactions <sup>†</sup>	300 reactions <sup>†</sup>
2x LAMP Master Mix	640 μL x 2 †12.5 μL/1 reaction	640 μL x 6

## III Storage

Store at -20°C in dark place.

## IV Protocol

1. Thaw the reagents completely at room temperature. Vortex the reagents to dissolve uniformly and then spin down and leave them to stand on ice.
2. Prepare Pre Mix.

2x LAMP Master Mix	12.5 μL
10x LAMP Primer Mix <sup>††</sup>	2.5 μL
Template	~ 5.0 μL
Total	up to 25.0 μL

<sup>††</sup>10x LAMP Primer Mix :

16 μM FIP, 16 μM BIP, 2 μM F3 Primer, 2 μM B3 Primer, 8 μM Loop Primer F, 8 μM Loop Primer B, 10 mM Tris-HCl (pH 8.0), 1 mM DTT

3. Transfer the Pre Mix to device recommendation tube and mix by inverting and spin down. Set to the fluorescence detecting device and start measuring (DNA amplification monitoring and Annealing curve analysis or Melting curve analysis). Follow the setting method of the device's manual.
4. 60~68°C<sup>†††</sup> for 30 min → Annealing curve analysis or Melting curve analysis.  
<sup>†††</sup> The optimum reaction temperature is changed depending on primer design.
5. Set Fluorescent wavelength to ResoLight Dye or SYBR<sup>®</sup> Green I when using Real-time PCR device.